

Preparation

1 SUITABILITY EVALUATION
 · Allocate 1 hour

2 LAY OFF YOUR CONTACT LENSES**
 · Soft : 1 week
 · Hard / RGP : 2 weeks

*** Subject to change based on the attending surgeon*


3 PRE-OPERATIVE SPECIALIST EYE EXAMINATION
 · Allocate about 2 - 3 hours
 · Arrange for transportation as your vision will be blurry

4 SURGERY DAY
 · Allocate about 2 - 3 hours at our centre and arrange your transportation


5 SURGERY DAY
 · 1 day, 1 week, 1 month & 1 year

INCLUSIVE: ONE (1) SET OF MEDICATION
 · Antibiotic drops
 · Anti-inflammatory drops
 · Artificial tears
 · Sunvisor
 · Protective eye shield


Why Us




Over **26+ years** of track record since **1999**








FDA Approved LASIK & Visian ICL




Lifetime Commitment for No Blade Lasik



Our **B Corp Certification** reflects our dedication to social & environmental responsibility

from blur to clear™



[+] ICL

Implantable Contact Lens

Nothing is removed so no major healing or adaptation is required. 5 - 10 minutes is all that is needed per eye.

Procedure

- 3mm incision made on cornea edge
- ICL is inserted between the iris and your natural lens
- ICL is positioned in the right position

Key Benefits

- Biocompatibility**
Collagen based material accepted by our body
- Predictable Results**
Nothing removed = no adaptation. Enhancement <0.2%[^] of cases
- High Range of Correction**
Up to -20.0 diopters (not limited by cornea thickness like laser surgeries)
- Permanent yet Reversible**
Designed to be permanent yet removable/replaceable if needed
- Small Wound Size**
- Excellent Night Vision**
- Safe**
- Quick Recovery**
- UV Protection**
- Flapless**

Technology by:

[-] Personalised LASIK

No Blade LASIK

A 100% bladeless, all-laser procedure tailored to your unique eye profile. It not only corrects short-sightedness and astigmatism but also enhances vision quality by reducing glare, halos, and other imperfections.

Procedure

- Individualized Vision Profile**
iDesign 2.0 Wavefront technology obtains a 3D map of each eye for a customized treatment of the visual imperfections
- Precise Corneal Flap**
A femtosecond laser is used to create a corneal flap
- Custom Laser Vision Correction**
Advanced CustomVue treatment is used to correct your vision

Key Benefits

- 25x Precision**
The iDesign 2.0 allows for 25x increased precision
- 5X Resolution**
Customized Treatment with 5X Resolution
iDesign 2.0 3D wavefront map allows the laser to treat imperfections in your eye with 5X the resolution
- Potentially Better than Glasses**
With the personalised treatment, you can potential see better vision, reduce poor night vision, glare & haloes
- No-Blade LASIK**
Enhance healing & reduce flap complication

Technology by:

[-] CLEAR PRO

The New Generation of Lenticular Extraction

A 2nd generation lenticular surgery that is **less invasive** and **minimizing your dry eyes**, leaving you to enjoy the more active lifestyle that you desire!

Procedure

- A thin layer of corneal tissue (lenticule) and a small incision is created with a low energy femtosecond laser
- The lenticule is removed with minimal stress to the cornea via the incision
- The resulting change of cornea shape corrects your problem and improves your vision

Key Benefits

- Flapless**
Enjoy an active lifestyle as there are no worries about flap complications or dislocations with CLEAR Pro!
- Less Dry Eye**
Being Flapless, CLEAR Pro helps preserve cornea biomechanical stability, while minimizing postoperative dry eye problems!
- Quick Recovery**
Low levels of energy pulses leads to a quicker wound recovery
- Small Wound Size**
With only a 2mm incision for CLEAR Pro, enjoy a faster wound recovery and lower infection risks!
- Safe**
2nd Generation Lenticule Extraction with AI assisted processes and eye tracking including cyclotorsion, via LIVE OCT monitoring for an accurate and safe procedure for best results

Technology by:

Find your **PERFECT MATCH**

	ICL	Personalised iLasik	CLEAR PRO
1 Addition / Removal#	Add [+]	Remove [-]	Remove [-]
2 Personalised	✓	✓✓	✗
3 Flapless	✓	✗	✓
4 Incision	Small	Big	Small
5 Suitable for dry eye	✓✓	✗	✓
6 Reversible	✓	✗	✗
7 Suitable for	Extended power range ◆◆ Dry eyes ◆◆ Prefer reversible treatment	Prefer customized treatment	Extreme sports/ Active lifestyle ◆◆ Dry eyes ◆◆ Worry flap complications
8 Lifetime Enhancement	✗	✓	✗
Pricing (per eye)	RM 7,998	RM 6,998	RM 6,998

- Standard Lasik: RM3,698 per eye
- Toric ICL (with astigmatism correction): RM9,998 per eye

#The surgery involves either the addition of a lens or the removal of corneal tissue.



Up to
24 Months
INSTALLMENT